

NEB EXPRESSIONS

A scientific update from New England Biolabs

Spring Edition 2011

Understanding Variability in
DNA Amplification Reactions *page 3*

OneTaq™ DNA Polymerase *page 6*

EpiMark™ Nucleosome
Assembly Kit *page 7*

HF Restriction Enzymes
and the NEBuffer advantage *pages 9*



Table of Contents

Special Offers

OneTaq™ iPad® Contest page 2

Feature Article

Understanding Variability
in Amplification Reactions page 3

Products

OneTaq™
DNA Polymerase..... pages 5-6

EpiMark™ Nucleosome
Assembly Kit page 7

HF Restriction Enzymes:
The NEBuffer advantage page 8-9

NEB's DNA Ladders..... page 10

NEBNext® Reagents for
Small RNA..... page 11

The 2011-12 NEB Catalog .. page 12

Upcoming Tradeshows

Visit the NEB booth at
the following meetings:

- ESHG, European Human Genetics Conference 2011
May 28 - 31, 2011
Amsterdam RAI, The Netherlands,
- EMBO Conference Series:
Chromatin and Epigenetics
June 1 - 5, 2011
EMBL Heidelberg, Germany
- 36th FEBS Congress (Federation of European Biochemical Societies)
June 25 - 30, 2011
Torino (Turin), Italy

 Follow NEB on Twitter

A Letter from NEB

Dear Researcher,

For over 30 years, New England Biolabs, Inc. has been at the forefront of the isolation, characterization and cloning of polymerases. Our experience and understanding of these integral reagents allows us to develop new polymerases that meet the needs of our customers, as well as optimize their use in PCR. The importance of understanding the interplay between variables in a PCR experiment is explored in the feature article. We also introduce our PCR Test Panel, a systematic approach designed to better understand the contribution of the numerous reagents and conditions in an amplification reaction.

The PCR Test Panel has helped us develop a new polymerase that is ideal for routine or problematic amplicons. OneTaq™ DNA Polymerase offers a robust solution to the vast majority of endpoint PCR needs and shows exceptional performance against other commercially available polymerases, especially when amplifying difficult or GC-rich templates.

Wishing you continued success in your research,

New England Biolabs



Spring at New England Biolabs' campus.

Try OneTaq™ for a Chance to Win an iPad®

Try our new OneTaq™ DNA Polymerase! Send us your feedback, including data, and receive a free NEB lab timer and a chance to win one of six iPads®*.

For contest rules, visit
www.neb.com/OneTaqContest

This offer is valid globally, while supplies last. One submission per person.

Submit early to increase your chances of winning. iPad winners will be drawn weekly on March 14, March 21, March 28 and April 4, and then monthly on May 2 and June 6, 2011.

iPad® is a registered trademark of Apple Computer, Inc.



Understanding Variability in DNA Amplification Reactions

The polymerase chain reaction (PCR, 1) is arguably the most common technique of molecular biology. As such, it has significantly impacted research and development in fields such as biochemistry, medicine, bioengineering and beyond. Success or failure in PCR is influenced by a myriad of factors including primer design, cycling conditions, and the quality and concentration of reaction substrates and solutions. By understanding the interplay of these variables, PCR-based tools and techniques are bolstered, and difficult amplifications become routine.

Nicole M. Nichols, Ph.D., New England Biolabs

The range of variables that impact amplification reactions is masked by the high success rate of typical PCR experiments. It is only after unsuccessful amplification attempts that these variables become evident. Ideally, various components can be altered to achieve success and even to favor a desired outcome (e.g., specificity over yield, sensitivity over specificity, etc.). Practical issues, such as very high (or low) template GC content, the presence of inhibitors, limitations of primers, source materials or time can limit the ability either to follow optimized PCR guidelines or to systematically evaluate a sufficient number of variables to ensure success.

The PCR Test Panel:

With the long-standing goal of enabling the research of our customers and our own scientists, NEB continues to devote resources both to basic DNA polymerase research and to the development of new amplification products. As part of these efforts, we have created a quantitative, microfluidic-based, PCR Test Panel. The goals of the PCR Test Panel are to:

1. Systematically manipulate the numerous variables present in an amplification reaction.
2. Understand the contributions of each variable to the desired signal-to-noise outcome.
3. Facilitate the development of enhanced tools (e.g., polymerases, buffers, etc.) that will broaden the definition of “routine PCR” to encompass situations that are currently challenging.

The PCR Test Panel is enabled by a microfluidic, agarose gel-based mimic (LabChip® GX platform, Caliper Life Sciences, Hopkinton, MA) that allows rapid end-point quantitation of amplification reactions in a 96-well plate format. A supporting, internal database has been designed to link the results from the microfluidic analysis (e.g., yield and purity of each expected product) to the detailed contents of each well (e.g., identity and concentration of polymerase, template, buffer, additives). Cycling conditions, thermocycler ID and other relevant details are also tracked.

By varying large numbers of conditions and quantitating both specificity and yield of a PCR prod-

uct, the contribution of each reaction component can be systematically evaluated. Some of the variations in template and primer components that are assessed by the PCR Test Panel are described in Table 1. Importantly, this system is easily adapted and can be scaled to assess new components and conditions as they arise.

In practice, the PCR Test Panel involves a typical workflow from experimental planning to data analysis, with crucial steps of database integration. The workflow is comprised of the following steps:

1. Design PCR experiment and capture relevant details in tracking database.
2. Run PCR experiments in triplicate (various module/experiment combinations to make 96-well plates).
3. Run 96-well plates through quantitative microfluidic analysis.
4. Export results from microfluidic platform into database to link quantitation with experimental details.


Experiments are performed in triplicate and relevant information about the contents of each well

is captured in the database. The data tracking system communicates the expected product sizes for each well/reaction to the microfluidic system. The output of microfluidic analysis is an electropherogram, which for familiarity is also represented as a typical gel-based image (Figure 1). The instrument software identifies the expected peaks and quantitates the yield and percent purity of each (described in more detail below). Finally, results from the analysis are imported into the database and connected with well contents and additional experimental information, forming the foundation of an expanding data warehouse of PCR results. Simple data analysis can be accomplished within the database, and IGOR Pro (Wavemetrics, Inc., Portland, OR) has been used for more complex analysis and graphical display.

Advances in Research and Development as a result of the PCR Test Panel:

The PCR Test Panel has furthered our understanding of many aspects of amplification reactions, one example being the best use for Lambda PCR tests. The functional tests of many commercially available polymerases employ Lambda as a substrate. However, the Lambda amplification

Figure 1. Electropherogram and interpreted gel image of primary PCR Test Panel data.



The electropherogram (left) shows overlaid results from two samples (G1, G4) selected from the interpreted gel (right). Expected peaks (EP) are identified by the software from information provided by the database. For product quantitation and size assessment, a DNA ladder is run after every 12 samples and upper and lower markers (UM, LM) are mixed 1:1 with each sample immediately prior to analysis.

Table 1: Examples of Test Panel Variables

| VARIABLE | EXAMPLES |
|-----------------|--|
| Template Source | Genomic DNA from simple (e.g., Lambda) and complex (e.g., plant, mouse, human) organisms |
| Length | From <100 bases to over 10 kb |
| GC/AT Content | From ~80% GC to ~ 80% AT |
| Repeats | Homopolymers, CpG islands, triplet repeats, etc. |
| Primer design | Typical and suboptimal primer pairs |

reactions that were part of the PCR Test Panel were so consistently robust that they did not serve as good indicators for how a polymerase would perform with real world, complex templates, even when matched for GC content. Interestingly, at typical template concentrations for Lambda, an increase in enzyme concentration resulted in a linear increase in the yield of amplified product (akin to enzyme titer experiments). A similar, linear response was not observed at working concentrations for more complex genomic targets \leq 1kb. Instead, high variation in yield was detected for individual amplicons in response to increased enzyme concentration in the reaction. These results demonstrate that Lambda can be used as an appropriate template for quality control assays, whereas templates relevant to daily laboratory experiments are better suited for optimizing product use recommendations.

Recently, the PCR Test Panel was used in the development of the *OneTaq* line of amplification

products. Extensive testing was used not only to guide formulations and eventual usage recommendations, but also to assess performance relative to a broad field of competitors. These comparisons demonstrated that *OneTaq* polymerases and buffers performed well in “routine” amplifications, but unlike many of their peers, continued to offer robust performance with more difficult amplicons, such as those with high GC or AT content. Figure 2 shows the gel-based image of triplicate reactions with *OneTaq* Hot Start DNA Polymerase (with and without the High GC enhancer) and six other hot start polymerases (plus competitor’s GC enhancers where provided) on a single high GC amplicon (68%). In addition, the quantitative primary output of microfluidic analysis (yield as ng/ μ l and % product purity as a function of the total signal in each reaction lane) allows a more condensed view of the data, which is more compatible with the scale of the PCR Test Panel. In Figure 3 (next page), the triplicate reactions shown in Figure 2 were condensed into a single column (boxed area in Figure 3) and are shown as part of a series of reactions with high GC, human genomic amplicons (ranging from 66%-80% GC content). One routine amplicon (55% GC) was also included as a control. Amplifications were all performed according to each manufacturer’s specific recommendations. Percent product purity is represented as circles in Panel A. Product yield is expressed as a bar chart in Panel B. Reactions were set up in the absence (solid bars) and presence (striped bars) of GC enhancers, where provided by the manufacturer. Percent purity scores account for contributions from primer-dimer formation and any other non-specific products

formed by the reaction. Low purity scores can arise from a prominent secondary product, from a smear of non-specific products, or from significant primer dimer interference (compare boxed area in Figure 3 to Figure 2 which shows the gel-based representation of the same data).

Although it would be convenient to suggest a single solution to all PCR difficulties, the reality is more complex and nuanced. For example, considering only GC content, we observed that the percent GC content of a template was not always an adequate predictor of reaction difficulty, nor could it alone define the need for a GC-specific buffer or enhancer. One striking example can be seen in Figure 3, where a 73% GC amplicon could only be robustly amplified using *OneTaq* in GC Buffer with the High GC Enhancer, whereas an 80% GC amplicon did not require the enhancer. Trends that emerge from studies of this scale help inform product guidelines for a broad range of product applications. For *OneTaq* polymerase, these include guidelines for increasing enzyme concentration when amplifying products over 3 kb and the use of a 68°C extension temperature. More information on *OneTaq* can be found on page 6 and at www.neb.com/OneTaq.


Conclusion:

Testing new products under a wide variety of conditions on a multitude of primer/template sets provides a comprehensive view of our products (and others in the market). Although it is interesting to study all the variables that affect PCR, we also realize that successful amplification of a desired target is what is important to our customers. The scope and scale of the PCR Test Panel allows us to see past the natural variability produced by the rugged PCR landscape where a single amplicon or reaction condition can be found to prove nearly any point desired. Data at this scale serves as a vantage point to build upon the strengths of various polymerases, buffers and protocols to expand the useful range of PCR conditions and the definition of what is routine. It allows our researchers at NEB to develop solutions to challenging amplifications and support your future polymerase-based applications.

References:


1. Saiki, et al., (1988) *Science* 239, 487-491
LabChip® is a registered trademark of Caliper Life Science.

Figure 2. Gel view of triplicate reactions on a single 68% GC amplicon for *OneTaq* Hot Start DNA Polymerase and hot start competitors (lanes #1 to #6).



A familiar gel output provides visual information about product yield and purity. Data from these triplicate reactions were averaged to create a single panel (*) of the graphs shown in Figure 3. Amplifications were all performed according to each manufacturer’s specific recommendations.

Figure 3: Examples of data from the PCR test panel: Comparison of OneTaq Hot Start DNA Polymerase to other commercially available hot start polymerases on high GC amplicons.



Reactions containing high GC human genomic DNA templates were set up at room temperature. PCR experiments included 30 cycles. Purity (A) and Yield (B) were calculated via microfluidic analysis from triplicate reactions. OneTaq polymerase was used in the absence (brown, solid) or presence (brown, striped) of High GC Enhancer. Competitor polymerases were cycled according to manufacturer's recommendations and included GC enhancers when supplied (striped bars).

Introducing

OneTaq™ DNA Polymerase

An optimized blend of *Taq* and Deep Vent_r DNA polymerases, OneTaq™ and OneTaq Hot Start DNA Polymerases offer robust amplification across a wide range of templates. The 3'–5' exonuclease activity of Deep Vent DNA Polymerase increases the fidelity and robustness of *Taq*, and the hot start formulation combines convenience with decreased interference from primer-dimers and secondary products. Available in convenient product formats, including master mixes, OneTaq shows exceptional performance against other commercially available polymerases, especially when amplifying difficult or GC-rich templates.

Advantages

- Exceptional performance in endpoint PCR across a wide range of templates
- Robust yields with minimal optimization
- Convenient product formats (stand-alone enzyme, master mixes, and Quick-Load® formats)
- Hot start version allows room temperature reaction setup and does not require a separate activation step
- Compatible with standard *Taq* protocols

See page 6 for more information




New Products

OneTaq: Robust amplification across a wide range of templates


OneTaq DNA Polymerase is supplied with two 5X buffers (Standard and GC), as well as a High GC Enhancer solution. For most routine and/or AT-rich amplicons or complex amplicons with up to ~65% GC content, OneTaq Standard Reaction Buffer provides robust amplification. For GC-rich amplicons, the OneTaq GC Reaction Buffer can improve both performance and yield. For particularly high GC amplicons (>65%) or difficult amplicons, the OneTaq High GC Enhancer can be added to reactions containing OneTaq GC Buffer. These formulations ensure maximum performance for routine, AT- or GC-rich amplicons.

Comparison of OneTaq DNA Polymerase to other commercially available polymerases (non-hot start).



Amplification of a selection of high GC human genomic DNA templates demonstrates OneTaq performance. PCR experiments included 30 amplification cycles. Purity (A) and Yield (B) were calculated via microfluidic analysis from triplicate reactions. Competitor polymerases were cycled according to manufacturer's recommendations.

Achieve robust amplification for routine, AT- and GC-rich templates with OneTaq



Amplification of a selection of sequences with varying AT and GC content from human and *C. elegans* genomic DNA using OneTaq DNA Polymerase. GC content is indicated above gel. Marker M is the 1 kb DNA Ladder (NEB #N3232).


Ordering Information

| PRODUCT | NEB # | SIZE |
|--|------------|-------------------------------|
| OneTaq™ DNA Polymerase | M0480S/L/X | 200/1,000/5,000 units |
| OneTaq™ 2X Master Mix with Standard Buffer | M0482S/L | 100/500 reactions (50 μl vol) |
| OneTaq™ 2X Master Mix with GC Buffer | M0483S/L | 100/500 reactions (50 μl vol) |
| OneTaq™ Quick-Load® 2X Master Mix with Standard Buffer | M0486S/L | 100/500 reactions (50 μl vol) |
| OneTaq™ Quick-Load® 2X Master Mix with GC Buffer | M0487S/L | 100/500 reactions (50 μl vol) |
| OneTaq™ Hot Start DNA Polymerase | M0481S/L/X | 200/1,000/5,000 units |
| OneTaq™ Hot Start 2X Master Mix with Standard Buffer | M0484S/L | 100/500 reactions (50 μl vol) |
| OneTaq™ Hot Start 2X Master Mix with GC Buffer | M0485S/L | 100/500 reactions (50 μl vol) |
| OneTaq™ Hot Start Quick-Load® 2X Master Mix with Standard Buffer | M0488S/L | 100/500 reactions (50 μl vol) |
| OneTaq™ Hot Start Quick-Load® 2X Master Mix with GC Buffer | M0489S/L | 100/500 reactions (50 μl vol) |

EpiMark™ Nucleosome Assembly Kit

NEB's line of EpiMark™ validated reagents for epigenetics studies now includes the EpiMark Nucleosome Assembly Kit. This kit contains the components necessary to form an unmodified recombinant human nucleosome using your own target DNA or the supplied control DNA. The protocol requires the mixing of already formed and purified recombinant human Histone H2A/H2B Dimer and Histone H3.1/H4 Tetramer in the presence of DNA at high salt, followed by dialysis down to low salt to make nucleosomes. One tetramer associates with two dimers to form the histone octamer on the DNA, generating a nucleosome. A method for assaying nucleosome formation by gel shift assay is also provided. These nucleosomes may serve as a better substrate for enzymes that are inactive on the DNA or one of the core histones alone. Each described reaction creates nucleosomes from ~50 pmol of a 208 bp DNA and may be scaled depending on the experiment.

Gel shift assay to visualize nucleosome assembly



Samples from nucleosome assembly reactions were run on 6% polyacrylamide gel in 0.5X TBE.

Lane 1: Low Molecular Weight DNA Ladder (NEB #N3233)

Lane 2: Nucleosome Control DNA

Lane 3: 0.5:1 ratio of Octamer* to DNA

Lane 4: 1:1 ratio of Octamer* to DNA

*Octamer = 2:1 mix of Histone H2A/H2B Dimer and Histone H3.1/H4 Tetramer.

The EpiMark Nucleosome Assembly Kit Includes:

- Histone H2A/H2B Dimer
- Histone H3.1/H4 Tetramer
- Control DNA

Ordering Information

| PRODUCT | NEB # | SIZE |
|---|--------|--------------|
| EpiMark™ Nucleosome Assembly Kit | E5350S | 20 reactions |
| COMPANION PRODUCTS | | |
| Histone H3.1/H4 Tetramer Human, Recombinant | M2509S | 1 nmol |
| Histone H2A/H2B Dimer Human, Recombinant | M2508S | 2 nmol |
| Nucleosome Control DNA | N1202S | 0.2 nmol |


Advantages

- Highly pure, recombinant system
- Pre-formed histone dimer and tetramer complexes simplify octamer formation
- Components stable for one year
- Dilution protocol only requires a few hours for assembly

Applications

- ChIP assay
- HAT assay
- Enzyme modification assays (i.e., methylation studies)

Purity of NEB Histones



Visit www.epimark.com to learn about additional products for epigenetics research available from NEB, including our broad line of recombinant human histones.



New Products


High Fidelity (HF) Restriction Enzymes


New England Biolabs provides customers with high quality tools for a wide range of molecular biology applications. As part of our ongoing commitment to the study and improvement of restriction enzymes, we are pleased to offer a line of restriction enzymes that have been engineered for maximum performance, convenience and flexibility.

EcoRI, BamHI, NotI, KpnI, HindIII, ...these restriction enzymes are some of the most commonly used workhorses of today's molecular biologists. However, under certain conditions these enzymes are reported to exhibit star activity. Star activity, or relaxed specificity, is an intrinsic property of restriction enzymes that leads to unwanted cleavage. In addition to its canonical consensus sequence "GAATTC" EcoRI for example may cut DNA at noncanonical sites ("NAATTC"), leading to non-predictable banding pattern and misinterpretation of experimental data. Different techniques such as cloning, genotyping, mutational analysis, mapping, probe preparation, sequencing and methylation detection employ a wide range of reaction conditions and often require the use of enzymes under suboptimal conditions.


As part of our ongoing commitment to the study and improvement of restriction enzymes, we are pleased to offer a superior line of High Fidelity (HFTM) restriction enzymes, that are compatible with our existing "wildtype" restriction enzymes:

EcoRI-HFTM is one of 20+ currently available optimized High Fidelity restriction enzymes. It has been engineered by NEB's scientists by exchanging functional amino acid residues . However, EcoRI-HFTM has the same specificity as the original established enzyme but with 64-fold reduced star activity .

All engineered High Fidelity (HFTM) enzymes offer the highest level of buffer compatibility, i.e. they work optimally in NEBuffer 4 . This will simplify double digest reactions, since more than 160 of NEB's restriction enzymes are supplied with NEBuffer 4 - the largest selection of restriction enzymes supplied in a single optimal buffer!






HFTM restriction enzymes have been engineered for robustness and flexibility to accelerate your research: all HFTM enzymes are also Time-SaverTM Qualified  and allow for an optional 5 min digest as well as overnight incubations.

To distinguish these engineered enzymes, the letters -HFTM have been added to the restriction enzyme name and they are packaged in unique purple-capped tubes.

NEB offers the new High Fidelity (HFTM) restriction enzymes at the same price as the original enzymes  - We recommend using HFTM restriction enzymes whenever possible in your DNA digests!

High Fidelity (HFTM) Restriction Enzymes -

Advantages:

-  Engineered for performance
-  Dramatically reduced Star Activity - Up to 64,000-fold reduced star activity compared to wildtype enzyme
-  Highest Buffer Compatibility - NEBuffer 4 is the optimal reaction buffer for over 160 restriction enzymes, including all HFTM enzymes
-  Speed - HFTM enzymes are Time-SaverTM qualified for 5 min digests
-  Value - HFTM enzymes are offered at the same fair price as the original enzymes

Use HFTM Restriction Enzymes in all applications for optimal results and convenience!

| Enzyme | optimal NEBuffer [†] | HF-Factor: reduced Star Activity compared to Wildtype* |
|--------------------------|-------------------------------|--|
| AgeI-HF TM | 4 + BSA | ≥ 8x |
| BamHI-HF TM | 4 | ≥ 125x |
| BsaI-HF TM | 4 + BSA | ≥ 250x |
| DraIII-HF TM | 4 + BSA | ≥ 1,000x |
| EagI-HF TM | 4 | 2x |
| EcoRI-HF TM | 4 | 64x |
| EcoRV-HF TM | 4 | ≥ 64x |
| HindIII-HF TM | 4 | ≥ 2,000x |


| Enzyme | optimal NEBuffer [†] | HF-Factor: reduced Star Activity compared to Wildtype* |
|------------------------|-------------------------------|--|
| KpnI-HF TM | 4 | ≥ 62,500x |
| MfeI-HF TM | 4 | ≥ 16x |
| NcoI-HF TM | 4 | ≥ 530x |
| NheI-HF TM | 4 + BSA | ≥ 266x |
| NotI-HF TM | 4 + BSA | ≥ 16x |
| PstI-HF TM | 4 | 33x |
| PvuI-HF TM | 4 + BSA | ≥ 32x |
| PvuII-HF TM | 4 | 32x |

| Enzyme | optimal NEBuffer [†] | HF-Factor: reduced Star Activity compared to Wildtype* |
|-----------------------|-------------------------------|--|
| SacI-HF TM | 4 + BSA | ≥ 266x |
| SalI-HF TM | 4 | ≥ 8,000x |
| SbfI-HF TM | 4 | 32x |
| Scal-HF TM | 4 | 62x |
| SphI-HF TM | 4 | 250x |
| SspI-HF TM | 4 | 16x |
| StyI-HF TM | 4 | 125x |

[†] Wild type enzymes were tested in supplied buffer for comparisons.

* Ref.: Wei, H. et al (2008) Nucleic Acids Reseach 36, e50.


Dramatically reduced star activity with HF restriction enzymes



PvuII shows 100% activity in all four NEBuffers, but exhibits intrinsic star activity in buffers other than the recommended NEBuffer 2 or under extended incubation times (i.e. overnight). Star activity is not observed with PvuII-HF™, even in overnight digests! (20 µl reactions, 2 µl of enzyme and incubated at 37°C. Marker M = 1 kb DNA Ladder #N3232.

Updated free NEB Activity Chart Poster available!

Request your free updated “NEB Activity Chart for Restriction Enzymes” Poster (75 x 45 cm) at your local distributor



Unrivalled buffer compatibility of NEB Restriction Enzymes: Save your budget with NEB!

New England Biolabs provides the largest selection of restriction enzymes that exhibit 100% activity in a single buffer. In fact, 169 Restriction Enzymes are supplied with NEBuffer 4; over 190 Restriction Enzymes are 100% active in this buffer!

Of the Top20 most commonly used restriction enzymes NEB provides 18 enzymes with NEBuffer 4. This unrivalled buffer compatibility results in unmatched convenience in reaction set-up and great savings to your budget!

NEB



NEB provides 18 enzymes with a single buffer NEBuffer4. Two different buffer are supplied with the Top20.

Invitrogen/ Life Technologies:



Six different buffers supplied with the Top20.

Roche:



Six different buffers supplied with the Top20.

Promega:



Seven different buffers supplied with the Top20.

Thermo/Fermentas:



Eight different buffers supplied with the Top20.*

*FastDigest Enzymes excluded

Top20 Enzymes: BamHI-HF, BglII, DpnI, EcoRI-HF, EcoRV-HF, Hinc II, Hind III-HF, KpnI-HF, NcoI-HF, NheI-HF, NotI-HF, PstI-HF, PvuII-HF, SacI-HF, SalI-HF, SmaI, SpeI, SphI-HF, Xba I, XhoI

Save money through NEB's buffer compatibility

Example: Double Digest EcoRI/XbaI:


NEB:

EcoRI-HF/XbaI both exhibit 100% activity in NEBuffer 4 ^{NEB4}, i.e. each enzyme works at the regular 1-fold concentration to obtain 100% cleavage.

vs.

Thermo/Fermentas:

There is no buffer available in which both enzymes exhibit 100% activity! Therefore an expensive 2-fold or 4-fold excess of XbaI is needed to get the job done!




Screenshot Fermentas Webpage Feb. 2011

Featured Product

DNA Ladders

NEB offers a variety of DNA Ladders with sizes ranging from 25 bp to 10 kb for use in agarose gel electrophoresis.



For PCR applications involving small DNA fragments try the special designed PCR Marker. More DNA, RNA and Protein ladders as well as conventional markers can be found at www.neb.com.

Ordering Information

| PRODUCT | NEB # | SIZE |
|--|----------|---------------------|
| CONVENTIONAL LADDERS AND PCR MARKER | | |
| 1 kb DNA Ladder | N3232S/L | 200/1,000 gel lanes |
| 100 bp DNA Ladder | N3231S/L | 100/500 gel lanes |
| 2-Log DNA Ladder | N3200S/L | 100/500 gel lanes |
| 50 bp DNA Ladder | N3236S/L | 100/500 gel lanes |
| Low Molecular Weight DNA Ladder | N3233S/L | 100/500 gel lanes |
| PCR Marker | N3234S/L | 30/150 µg |
| Supercoiled DNA Ladder | N0472S | 100 gel lanes |
| QUICK-LOAD® DNA LADDERS | | |
| Quick-Load 2-Log DNA Ladder | N0469S/L | 125/375 gel lanes |
| Quick-Load 1 kb DNA Ladder | N0468S/L | 125/375 gel lanes |
| Quick-Load 100 bp DNA Ladder | N0467S/L | 125/375 gel lanes |
| Quick-Load 50 bp DNA Ladder | N0473S | 125 gel lanes |
| Quick-Load Low Molecular Weight DNA Ladder | N0474S | 125 gel lanes |
| Quick-Load PCR Marker | N0475S | 125 gel lanes |
| TRIDYE™ DNA LADDERS | | |
| TriDye 2-Log DNA Ladder | N3270S | 125 gel lanes |
| TriDye 1 kb DNA Ladder | N3272S | 125 gel lanes |
| TriDye 100 bp DNA Ladder | N3271S | 125 gel lanes |


Advantages

- Stable at room temperature
- Sharp uniform bands
- Easy-to-identify reference bands
- Can be used for sample quantitation
- No extra backbone DNA present

Convenient Formats

The **Quick-Load DNA Ladders** are pre-mixed, ready-to-load molecular weight markers containing bromophenol blue as a tracking Dye.







The **TriDye DNA Ladders** are pre-mixed, ready-to-load molecular weight markers containing 3 dyes which serve as visual aids to monitor the progress of migration during agarose gel electrophoresis.



TriDye™ Ladders: 1 kb, 50 bp, 100 bp and 2-Log DNA Ladders are available in ready-to-load format containing xylene cyanol FF, bromophenol blue and orange G.

Novel Protocol for Small RNA Sample Prep for Next Gen Sequencing

New England Biolabs has expanded its line of NEBNext reagents to include NEBNext Small RNA Sample Prep Set 1 for sequencing on the Illumina® platforms, and NEBNext Small RNA Sample Prep Set 3 for sequencing on Life's SOLiD platform. These new products offer a new and unique protocol, which results in higher yields and lower adaptor-dimer formation, as compared to other protocols. With oligos, enzymes and buffers included, both sets provide substantial cost savings, and the master mix format streamlines the workflow by reducing the number of vials and pipetting steps required.

| | NEBNext® Small RNA Sample Prep Set 1 Reagents Supplied (Illumina Compatible) | NEBNext® Small RNA Sample Prep Set 3 Reagents Supplied (SOLiD Compatible) |
|--|---|---|
| NEB # | E6120S/L | E6160S/L |
| Size | 10/50 reactions | 10/50 reactions |
| 3' Ligation  | <ul style="list-style-type: none"> • 3' Ligation Enzyme Mix • 3' Ligation Reaction Buffer (2X) • 3' SR Adaptor 1 | <ul style="list-style-type: none"> • 3' Ligation Enzyme Mix • 3' Ligation Reaction Buffer (2X) • 3' SR Adaptor 3 |
| Primer Hybridization  | <ul style="list-style-type: none"> • SR RT Primer 1 | <ul style="list-style-type: none"> • SR RT Primer 3 |
| 5' Ligation  | <ul style="list-style-type: none"> • 5' Ligation Enzyme Mix • 5' Ligation Reaction Buffer (10X) • 5' SR Adaptor 1 • Nuclease-Free Water | <ul style="list-style-type: none"> • 5' Ligation Enzyme Mix • 5' Ligation Reaction Buffer (10X) • 5' SR Adaptor 3 • Nuclease-Free Water |
| First Strand Synthesis  | <ul style="list-style-type: none"> • RNase Inhibitor, Murine • dNTPs | <ul style="list-style-type: none"> • RNase Inhibitor, Murine • dNTPs |
| PCR Amplification  | <ul style="list-style-type: none"> • LongAmp® Taq 2X Master Mix • SR Primer F1 • SR Primer R1 | <ul style="list-style-type: none"> • LongAmp® Taq 2X Master Mix • SR Primer F3 • SR Primer R3 |
| Size Selection  | <ul style="list-style-type: none"> • Gel Loading Dye, Orange (6X) • Quick-Load® Low Molecular Weight DNA Ladder • DNA Gel Elution Buffer (1X) • Linear Acrylamide • TE | <ul style="list-style-type: none"> • Gel Loading Dye, Orange (6X) • Quick-Load® Low Molecular Weight DNA Marker • DNA Gel Elution Buffer (1X) • Linear Acrylamide • TE |


| | | |
|------------|------------------|---------------------|
| KEY | RNA = Red | cDNA = Blue |
| | Adaptor = Green | DNA Product = Black |
| | RT Primer = Aqua | |
| | | |

For the complete listing of NEBNext reagents available, please visit www.nebnext.com

Advantages of Small RNA Sets

- **Convenient formats** – Enzymes, buffers and nucleotides are included at the appropriate concentrations and in appropriate volumes, and are available in set or master mix format
- **Novel Protocol for small RNA** – Higher yield and substantially reduced adaptor dimer formation suitable for methylated and unmethylated small RNAs
- **Functional Validation** – Each set is functionally validated by preparation of a library from a standard reference RNA, followed by Illumina or SOLiD sequencing
- **Stringent Quality Controls** – Additional QCs ensure maximum quality and purity
- **Value Pricing**

NEBNext Small RNA sample prep increases library yield and substantially reduces adaptor-dimer contamination



Standard and NEBNext protocols for small RNA Library construction were carried out using 5 µg of rat testis total RNA. Note: miRNAs as well as other small RNA species (including piRNAs that are very abundant in rat testis), are efficiently captured by the NEBNext protocol.

Now Available: The 2011–12 NEB Catalog and Technical Reference

The NEB Catalog & Technical Reference contains over 100 new products in areas such as DNA cloning, PCR, epigenetics, RNA analysis, sample prep for next generation sequencing, protein expression and cellular analysis. In addition, our popular technical reference section includes up-to-date technical charts, protocols and troubleshooting tips to aid experimental design.



Visit www.neb.com/litrequest
or contact your local distributor
and order your copy today!

Catalog highlights:

- An expanded line of High Fidelity (HF) Restriction Enzymes engineered for reduced star activity and convenience
- A comprehensive offering of PCR reagents
- A new section highlighting NEBNext reagents for sample preparation of DNA and RNA for next generation sequencing
- The latest innovation in competent cells for protein expression
- A new epigenetics section showcasing novel methods for the identification and quantitation of 5-mC and 5-hmC within a DNA locus
- Tools to enable RNA research
- An extensive selection of markers and ladders for DNA, RNA and protein analysis
- A powerful protein labeling technology unique to NEB
- A range of expression systems, including a novel kit for cell-free expression
- Tools for glycobiology



DNA CLONING
DNA AMPLIFICATION & PCR
EPIGENETICS
RNA ANALYSIS
SAMPLE PREP FOR NEXT GEN SEQUENCING
PROTEIN EXPRESSION & ANALYSIS
CELLULAR ANALYSIS

Your local NEB distributor:

BeNeLux:
BIOKÉ
Plesmanlaan 1d
2333 BZ Leiden
The Netherlands
Tel: +31(0)71 568 1000
Tel: Belgium: 0800 71640
Fax: +31(0)71 568 1010
info@bioke.com
www.bioke.com

ITALY:
EUROCLONE S.P.A.
Via Figino 20/22
20016 Pero (Milan)
Free call: 800-315911
Tel: (02) 381951
Fax: (02) 38101465
info@euroclone.it
www.euroclonergroup.it/celbio/

CZECH REPUBLIC:
BIOTECH a.s.
Sluzeb 4
108 52 Praha 10
Tel: (Toll Free) 0800 124683
Fax: 02 72701742
biotech@biotech.cz
www.biotech.cz

JORDAN:
BURJOUAN
Est. for Medical & Lab Supplies
Jreisat Center 2nd floor
P.O.Box 2425
11910 Amman
ahmad-m@burjouan.com
www.burjouan.com

DENMARK:
MEDINOVA SCIENTIFIC A/S
Produktionsvej 26
2600 Glostrup
Tel: (39) 56 20 00
Fax: (39) 56 19 42
medinova@medinova.dk
www.medinova.dk

POLAND:
LENCOMM TRADE
INTERNATIONAL
ul.Wolczynska 133
01-919 Warszawa
tel:(+48 22) 835-62-78
fax:(+48 22) 835-03-54
office@lencomm.pl
www.lencomm.pl

FINLAND:
FISHER SCIENTIFIC
Ratastie 2
01620 Vantaa
Tel: +358 9 329 100
Fax: +358 9 3291 0500
finnzymes@thermofisher.com
www.fishersci.fi

Lab-JOT Tomasz Dzieczewski
ul. Nadarzyńska 45a, Kanie,
05-805 Orrebusy
stacjonarny: (+48 22) 758 59 07
komórkowy: (+48 606) 338 879
fax: (+48 22) 758 59 07
labjot@gmail.com
www.labjot.com

FRANCE:
OZYME
BP 268
78053 Saint Quentin Yvelines CEDEX
Tel: (1) 34 60 24 24
Fax: (1) 34 60 92 12
info@ozyme.fr
www.ozyme.fr

SWEDEN:
IN VITRO SWEDEN AB
Box 21160
S-100 31 Stockholm
Tel: (08) 30 60 10
Fax: (08) 30 60 15
iv@invitro.se
www.invitro.se

GREECE:
BIOLINE SCIENTIFIC
1 Meg. Alexandrou Str.
10437 Athens
Tel: 210 5226547
Fax: 210 5244744
demagos@hol.gr

SPAIN/PORTUGAL:
IZASA S.A.
Aragón 90
08015 Barcelona
Tel: 902.20.30.70
Fax: 902.20.30.71
consultasbiotec@izasa.es
www.izasa.com

HUNGARY:
KVALITEX SCIENTIFIC
Technological Trading Ltd.
Pannónia u. 7.
1136 Budapest
Phone: (1) 340-4700
Fax: (1) 339-8274
info@kvalitex.hu
www.kvalitex.hu

SWITZERLAND:
BIOCONCEPT
Innovationszentrum
Nordwestschweiz
Paradiesrain 14
Postfach 427
CH 4123 Allschwil 1
Tel: (061) 486 80 80
Fax: (061) 486 80 00
info@bioconcept.ch
www.bioconcept.ch

ISRAEL:
ORNAT BIOCHEMICALS AND
LABORATORY EQUIPMENT LTD.
POB 2071
Rehovot 76120
Tel: +972-8-9477077
Fax: +972-8-9363934
ornatbio@ornat.co.il
www.ornat.co.il

TURKEY:
SACEM LIFE TECH
M. Fevzi Cakmak Cad. 8/4
06500 Bahcelievler Ankara
Tel: +90/312/223 52 72
Fax +90/312/215 48 80
sacem@sacem.com.tr
www.sacem.com.tr