



# Circulating DNA from plasma

## User manual

NucleoSpin<sup>®</sup> 96 DNA Plasma  
NucleoSpin<sup>®</sup> 96 DNA Plasma Core Kit

August 2016 / Rev. 01

## Table of contents

1	Components	4
1.1	Kit contents	4
1.2	Reagents, consumables, and equipment to be supplied by user	5
1.3	About this user manual	6
2	Product description	7
2.1	The basic principle	7
2.2	Kit specifications	7
2.3	Size and yield of DNA from plasma	8
2.4	Handling of sample material	8
2.5	Elution procedures	8
2.6	Stability of isolated DNA	8
3	Storage conditions and preparation of working solutions	9
4	Safety instructions	10
5	Protocol for the isolation of DNA from plasma	13
6	Appendix	17
6.1	Troubleshooting	17
6.2	Ordering information	18
6.3	Product use restriction /warranty	18

# 1 Components

## 1.1 Kit contents

<b>NucleoSpin® 96 DNA Plasma</b>		
<b>REF</b>	<b>1 x 96 preps 740873.1</b>	<b>4 x 96 preps 740873.4</b>
Activation Buffer PMA	75 mL	2 x 100 mL
Lysis Buffer PML	100 mL	4 x 100 mL
Binding Buffer PMB	2 x 250 mL	8 x 250 mL
Wash Buffer PMW1	100 mL	500 mL
Wash Buffer PMW2 (Concentrate) *	50 mL	2 x 100 mL
Elution Buffer PME **	13 mL	60 mL
Liquid Proteinase K	6 mL	4 x 6 mL
NucleoSpin® Plasma Binding Plate	1	4
MN Wash Plate	1	4
Square-well Block (96 wells)	2	8
Round-well Block with Cap Strips	1	4
User manual	1	1

\* For preparation of working solutions and storage conditions see section 3

\*\* Elution Buffer PME: 5 mM Tris/HCl, pH 8.5

<b>NucleoSpin® 96 DNA Plasma Core Kit</b>		
<b>REF</b>	<b>1 x 96 preps 740874.1</b>	<b>4 x 96 preps 740874.4</b>
Activation Buffer PMA	75 mL	2 x 100 mL
Lysis Buffer PML	100 mL	4 x 100 mL
Binding Buffer PMB	2 x 250 mL	8 x 250 mL
Wash Buffer PMW1	100 mL	500 mL
Wash Buffer PMW2 (Concentrate) *	50 mL	2 x 100 mL
Elution Buffer PME	13 mL	60 mL
Liquid Proteinase K	6 mL	4 x 6mL
NucleoSpin® Plasma Binding Plate	1	4
User manual	1	1

## 1.2 Reagents, consumables, and equipment to be supplied by user

### Reagents:

- 96–100 % ethanol

### Consumables:

- For lysis and adjusting of binding conditions: 1 x 96 Square-well Block for processing 500 µL samples, 2 x 96 Square-well Blocks for processing 1 mL samples, 4 x 96 Square-well Blocks for processing 2 mL samples. Alternatively, 24-Square-well Blocks can be used (see ordering information section 6.2).
- Disposable pipette tips

### Equipment:

- NucleoVac 96 Vacuum Manifold (see ordering information section 6.2)
- NucleoVac Vacuum Regulator (see ordering information section 6.2)
- Vacuum pump
- Heater-shaker or incubator oven
- Multi channel pipettes or large volume pipettes with appropriate tips
- Personal protection equipment (lab coat, gloves, goggles)

\* For preparation of working solutions and storage conditions see section 3

\*\* Elution Buffer PME: 5 mM Tris/HCl, pH 8.5

### **1.3 About this user manual**

It is strongly recommended for first time users to read the detailed protocol sections of the **NucleoSpin® 96 DNA Plasma** kit before using this product. Experienced users, however, may refer to the Protocol-at-a-glance instead. The Protocol-at-a-glance is designed to be used only as a supplemental tool for quick referencing while performing the purification procedure.

All technical literature is available online at [www.mn-net.com](http://www.mn-net.com).

## 2 Product description

### 2.1 The basic principle

The **NucleoSpin® 96 DNA Plasma** kit is designed for the efficient isolation of circulating DNA from human blood plasma. Fragmented DNA of 50 bp and larger can be purified with high efficiency. The **NucleoSpin® 96 DNA Plasma** kit can be used with manual and automated vacuum manifolds. The kit is fully automatable on many liquid handling robots.

The protocol follows state-of-the-art bind-wash-elute procedures: lysis is performed within 30 minutes with Proteinase K and lysis buffer. Afterwards, a binding buffer is added and the solution is applied onto the wells of the binding plate in several steps and DNA is bound to the silica membrane. Three washing steps efficiently remove contaminating substances, such as PCR inhibitors. Drying of silica is achieved by applying vacuum and pure DNA is finally eluted.

### 2.2 Kit specifications

- The **NucleoSpin® 96 DNA Plasma** kit is recommended for the isolation of circulating cell-free DNA from human plasma.
- The **NucleoSpin® 96 DNA Plasma** kit is designed for high recovery of fragmented DNA  $\geq$  50 bp.
- Up to 2 mL plasma can be used as sample material in a single well
- DNA yield strongly depends on the individual sample, but is typically in the range of 0.1 to 100 ng DNA per mL plasma.
- DNA is ready to use for downstream applications such as real-time PCR or NGS.
- The preparation time is approximately 90 min for up to 96 plasma samples.

#### Kit specifications at a glance

Parameter	NucleoSpin® 96 DNA Plasma
Technology	Silica-membrane technology
Format	96-well plates
Sample material	Human EDTA/Cell-Free DNA BCT® (Streck) plasma
Sample amount	0.5–2 mL per preparation
Typical yield	sample dependent
Elution volume	100 $\mu$ L (recovered volume about 70 $\mu$ L)
Preparation time	approx. 90 min / 96 preps

## 2.3 Size and yield of DNA from plasma

Usually, DNA concentrations in plasma are in a range of 0.1 ng up to several 100 ng DNA per mL of plasma. The amount of circulating DNA in plasma depends on health condition of the donor, sampling and handling of the blood, plasma preparation, DNA isolation method, and others.

A significant portion of the cell-free DNA in plasma originates from apoptotic cells. Therefore, a considerable percentage of this circulating DNA is known to be highly fragmented. However, the degree of fragmentation and the ratio of fragmented DNA to high molecular weight DNA depends on several parameters like origin of the DNA (e.g., fetal, tumor, microbial DNA), health condition of the blood donor, blood sampling procedure, and handling of the sample.

## 2.4 Handling of sample material

Circulating DNA yield and quality is largely influenced by blood sampling technique, handling, storage, and plasma preparation. It is highly recommended to perform these steps as uniform as possible in order to achieve highest reproducibility.

Plasma can be isolated according to the following recommendation:

### **Preparation of plasma from human EDTA blood or Streck Cell-Free DNA BCT®**

- 1 Centrifuge samples for 10 min at 2,000 x *g*.
- 2 Remove the plasma without disturbing sedimented cells and particles.
- 3 Clear plasma of residual cellular debris by means of centrifugation (10 min at 5000 x *g*).
- 4 If necessary, freeze plasma samples in fresh tubes. Upon thawing, check for precipitates and remove them with a final centrifugation step.

## 2.5 Elution procedures

The recommended standard elution procedure comprises two steps of 50 µL. This will result in about 70 µL eluate. The retained volume will contain very little amounts of DNA because the majority will be present in the eluted fraction.

## 2.6 Stability of isolated DNA

Due to the low DNA content in plasma and the resulting low total amount of isolated DNA, its fragmentation, and the absence of DNase inhibitors (the elution buffer does NOT contain EDTA), the eluates should be kept on ice for short term storage and frozen at -20 °C for long term storage.

### 3 Storage conditions and preparation of working solutions

**Attention:** Buffers PML, PMB and PMW1 contain guanidinium hydrochloride (chaotropic salt) which can form highly reactive compounds when combined with bleach (sodium hypochlorite). DO NOT add bleach or acidic solutions directly to the sample-preparation waste. Wear gloves and goggles!

Storage conditions:

- All kit components can be stored at room temperature (18–25 °C) and are stable for at least one year.
- If there is any precipitate present in the buffers, warm the buffer up to 25–37 °C to dissolve the precipitate before use.

Before starting any **NucleoSpin® 96 DNA Plasma** protocol, prepare the following:

- **Wash Buffer PMW2:** Ethanol has to be added to Wash Buffer PMW2 (Concentrate) according to the instructions on the label and in this user manual. Mark the label of the bottle to indicate that ethanol was added. All other kit components are ready to use.
- Prepare plasma sample according to section 2.4.
- Set heating block or incubator oven to 56 °C for lysis.
- Set up the **NucleoVac 96 Vacuum Manifold**.
- **Liquid Proteinase K** is ready to use. After first opening, store Liquid Proteinase K at -20 °C.
- When using multi well plates, samples have to be split into suitable aliquots.

#### NucleoSpin® 96 DNA Plasma

REF	1 x 96 preps 740873.1	4 x 96 preps 740873.4
Buffer PMW2 (Concentrate)	50 mL Add 200 mL ethanol	2 x 100 mL Add 400 mL ethanol to each bottle

#### NucleoSpin® 96 DNA Plasma Core Kit

REF	1 x 96 preps 740874.1	4 x 96 preps 40874.4
Buffer PMW2 (Concentrate)	50 mL Add 200 mL ethanol	2 x 100 mL Add 400 mL ethanol to each bottle



## 4 Safety instructions

The following components of the **NucleoSpin® 96 DNA Plasma** and **NucleoSpin® 96 DNA Plasma Core Kit** contain hazardous contents.

*Wear gloves and goggles and follow the safety instructions given in this section.*

### GHS classification

Only harmful features do not need to be labeled with H and P phrases up to 125 mL or 125 g.

*Mindergefährliche Eigenschaften müssen bis 125 mL oder 125 g nicht mit H- und P-Sätzen gekennzeichnet werden.*

Component	Hazard contents	GHS symbol	Hazard phrases	Precaution phrases
<i>Inhalt</i>	<i>Gefahrstoff</i>	<i>GHS-Symbol</i>	<i>H-Sätze</i>	<i>P-Sätze</i>
Activation Buffer PMA	Sodium hydroxide solution 0.5–2.0 % <i>Natriumhydroxid-Lösung 0,5–2,0 %</i> CAS 1310-73-2	 WARNING ACHTUNG	290, 315, 319	234, 264, 280, 302+352, 305+351+338, 332+313, 337+313, 390, 406
Lysis Buffer PML	Guanidine hydrochloride 50–66 % <i>Guanidinhydrochlorid 50–66 %</i> CAS 50-01-1	 WARNING ACHTUNG	302, 315, 319	264, 280, 301+312, 302+352, 305+351+338, 330, 332+313, 337+313
Binding Buffer PMB	Guanidine hydrochloride 24–36 % and ethanol 35–55 % <i>Guanidinhydrochlorid 24–36 % und Ethanol 35–55 %</i> CAS 50-01-1, 64-17-5	 WARNING ACHTUNG	226, 302	210, 233, 301+312, 330, 370+378, 403+235
Wash Buffer PMW1	Guanidine hydrochloride 36–50 % and 2-propanol 20–50 % <i>Guanidinhydrochlorid 36–50 % und 2-Propanol 20–50 %</i> CAS 50-01-1, 67-63-0	 WARNING ACHTUNG	226, 302, 319, 336	210, 233, 264, 280, 301+312, 305+351+338, 330, 337+313, 370+378, 403+235
Liquid Proteinase K	Proteinase K, liquid 1–3 % <i>Proteinase K, flüssig 1–3 %</i> CAS 39450-01-6	 WARNING ACHTUNG	317	261, 272, 280, 302+352, 333+313, 363

### Hazard phrases

H 226 Flammable liquid and vapor  
*Flüssigkeit und Dampf entzündbar.*

- H 290 May be corrosive to metals.  
*Kann gegenüber Metallen korrosiv sein.*
- H 302 Harmful if swallowed.  
*Gesundheitsschädlich bei Verschlucken.*
- H 315 Causes skin irritation.  
*Verursacht Hautreizungen.*
- H 317 May cause an allergic skin reaction.  
*Kann allergische Hautreaktionen verursachen.*
- H 319 Causes serious eye irritation.  
*Verursacht schwere Augenreizung*
- H 336 May cause drowsiness or dizziness.  
*Kann Schläfrigkeit und Benommenheit verursachen.*

### Precaution phrases

- P 210 Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking.  
*Von Hitze, heißen Oberflächen, Funken, offenen Flammen sowie anderen Zündquellenarten fernhalten. Nicht rauchen.*
- P 233 Keep container tightly closed.  
*Behälter dicht verschlossen halten.*
- P 234 Keep only in original container.  
*Nur im Originalbehälter aufbewahren.*
- P 261 Avoid breathing dust / fume / gas / mist / vapors / spray.  
*Einatmen von Staub / Rauch / Gas / Nebel / Dampf / Aerosol vermeiden.*
- P 264 Wash ... thoroughly after handling.  
*Nach Handhabung ... gründlich waschen.*
- P 272 Contaminated work clothing should not be allowed out of the workplace.  
*Kontaminierte Arbeitskleidung nicht außerhalb des Arbeitsplatzes tragen.*
- P 280 Wear protective gloves / protective clothing / eye protection / face protection.  
*Schutzhandschuhe / Schutzkleidung / Augenschutz / Gesichtsschutz tragen.*
- P 301+312 IF SWALLOWED: Call a POISON CENTER / doctor / ... / if you feel unwell.  
*BEI VERSCHLUCKEN: Bei Unwohlsein GIFTINFORMATIONSZENTRUM / Arzt / ... anrufen.*
- P 302+352 IF ON SKIN: Wash with plenty of water / ...  
*BEI BERÜHRUNG MIT DER HAUT: Mit viel Wasser / ... waschen.*
- P 305+351+338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
*BEI KONTAKT MIT DEN AUGEN: Einige Minuten lang behutsam mit Wasser ausspülen. Eventuell vorhandene Kontaktlinsen nach Möglichkeit entfernen. Weiter ausspülen.*
- P330 Rinse mouth.  
*Mund ausspülen.*
- P 332+313 If skin irritation occurs: Get medical advice/attention.  
*Bei Hautreizung: Ärztlichen Rat einholen/ärztliche Hilfe hinzuziehen.*
- P 333+313 If skin irritation or rash occurs: Get medical advice / attention.  
*Bei Hautreizung oder -ausschlag: Ärztlichen Rat einholen / ärztliche Hilfe hinzuziehen.*
- P 337+313 If eye irritation persists: Get medical advice/attention.  
*Bei anhaltender Augenreizung: Ärztlichen Rat einholen/ärztliche Hilfe hinzuziehen.*

- P 363 Wash contaminated clothing before reuse.  
*Kontaminierte Kleidung vor erneutem Tragen waschen.*
- P 370+378 In case of fire: Use ... to extinguish.  
*Bei Brand: ... zum Löschen verwenden.*
- P 390 Absorb spillage to prevent material damage.  
*Ausgetretene Mengen zur Vermeidung von Materialschäden aufnehmen.*
- P403+235 Store in a well-ventilated place. Keep cool.  
*An einem gut belüfteten Ort aufbewahren. Kühl halten.*
- P 406 Store in corrosive resistant / ... container with a resistant inner liner.  
*In korrosionsbeständigem / ... Behälter mit widerstandsfähiger Innenauskleidung aufbewahren.*



The symbol shown on labels refers to further safety information in this section.  
*Das auf Etiketten dargestellte Symbol weist auf weitere Sicherheitsinformationen dieses Kapitels hin.*

## 5 Protocol for the isolation of DNA from plasma

### Protocol-at-a-glance

<b>1 Lyse samples</b>	<b>25 µL Proteinase K</b>
	<b>1 mL plasma</b>
	<b>Mix</b>
	<b>400 µL Buffer PML</b>
	<b>Mix</b>
	<b>Incubate at 56 °C, 30 min</b>
	<i>Note: For Streck Cell-Free DNA BCT® incubate 60 min</i>
<b>Prepare the NucleoVac 96 Vacuum Manifold</b>	
	<b>400 µL Buffer PMA per well</b>
	<b>-0.4 bar*, 1 min</b>
<b>2 Adjust DNA binding conditions</b>	<b>2 mL Buffer PMB</b>
	<b>Mix</b>
<b>3 Bind DNA</b>	<b>Transfer lysates</b>
	<i>Note: Transfer lysates in aliquots of 1 mL</i>
	<b>-0.4 bar*, 2 min</b>
<b>4 Wash silica membrane</b>	<b>800 µL PMW1</b>
	<b>-0.4 bar*, 2 min</b>
	<b>1 mL PMW2</b>
	<b>-0.4 bar*, 2 min</b>
	<b>1 mL PMW2</b>
	<b>-0.4 bar*, 2 min</b>
<b>5 Dry silica membrane</b>	<b>-0.6 bar*, 10 min</b>

\* Reduction of atmospheric pressure

**6 Elute DNA**

**50 µL PME**

**Incubate 1 min at RT**

**-0.4 bar\*,  
30 s**


**50 µL PME**

**-0.6 bar\*,  
30 s**

**Setup of vacuum manifold:**

**Binding / Washing steps**

**Elution step**




**Step 4:**  
Place the NucleoSpin® Plasma Binding Plate on top of the manifold lid.

**Step 3:**  
Place the manifold lid on top of the manifold base.

**Step 2:**  
Place the MN Wash Plate in the manifold.

**Step 1:**  
Insert spacers 'MTP/MULTI-96 PLATE' and waste container in the manifold base.




**Step 4:**  
Place the NucleoSpin® Plasma Binding Plate on top of the manifold lid.


**Step 3:**  
Place the manifold lid on top of the manifold base.

**Step 2:**  
Place the Round-well Block in the manifold.

**Step 1:**  
Insert spacers 'ROUNDWELL BLOCK' in the manifold base.



**Final setup**



**Final setup**

\* Reduction of atmospheric pressure

## 5.1 Detailed protocol

The procedure below describes the isolation of cell-free DNA from 1 mL human plasma. Adjusting reagent volumes according to the table below allows for processing of plasma volumes from 0.5–2 mL.

Plasma volume [mL]	Liquid Proteinase K [μL]	Lysis Buffer PML [μL]	Binding Buffer PMB [mL]
0.5	12.5	200	1
1	25	400	2
2	50	800	4

### 1 Lyse sample

Add **25 μL Liquid Proteinase K** to each well of a 24-Square-well block. (Alternatively, samples can be split into aliquots of 500 μL and processed in 96 Square-well Blocks.)

Add **1 mL plasma** and **mix** by pipetting.

Add **400 μL Buffer PML** and **mix** by pipetting.

Incubate at **56 °C for 30 min** (For Streck Cell-Free DNA BCT<sup>®</sup>, incubate 60 min; ideally with shaking).

#### Prepare the NucleoVac 96 Vaccum manifold

Insert spacers 'MTP/MULTI-96 PLATE', the Waste Container and the Wash Plate into the manifold base. Place the manifold lid on top and then the **NucleoSpin<sup>®</sup> Plasma Binding Plate**.

While incubating the lysis, apply **400 μL Buffer PMA** to the wells of the NucleoSpin<sup>®</sup> Plasma Binding Plate. Incubate one minute, then apply vacuum of **-0.4 bar\*** for **1 min**.

### 2 Adjust binding conditions

Add **2 mL Buffer PMB** to each well and mix by pipetting.

### 3 Bind DNA

**Transfer** prepared lysates to the NucleoSpin<sup>®</sup> Plasma Binding Plate in aliquots of 1 mL.

Apply vacuum of **-0.4 bar\*** for **2 min**.

Lysates can be loaded continuously while they are passing the membrane. Remove, empty and replace Waste Container after the third loading step as well as after column loading is completed.

\* Reduction of atmospheric pressure

#### 4 Wash silica membrane

##### **1<sup>st</sup> wash**

Once all lysates have passed the membrane, add **800 µL Buffer PMW1** to each well. Incubate for **1 min**, and then apply vacuum of **-0.4 bar\*** for **2 min**.

##### **2<sup>nd</sup> wash**

Add 1 mL **Buffer PMW2**, and then apply vacuum of **-0.4 bar\*** for **2 min**.

##### **3<sup>rd</sup> wash**

Repeat 2<sup>nd</sup> wash step.

---

Remove and empty Waste Container. Remove Wash Plate.

---

#### 5 Dry silica membrane

Apply strongest possible vacuum of at least **-0.6 bar\*** for **10 min** to dry the silica membrane.

---

After drying, blot column outlets on tissue paper to remove residual ethanol.

---

Insert spacers "ROUNDWELL Block" and place a Round-well Block on top. Place the manifold lid on top and then the **NucleoSpin® Plasma Binding Plate**.

---

#### 6 Elute highly pure DNA

Add **50 µL Buffer PME** (first elution step) to the membrane. Incubate **1 min**.

Apply vacuum of **-0.4 bar\*** for **30 s**.

Add **50 µL Buffer PME** (second elution step) to the membrane.

Apply vacuum of **-0.6 bar\*** for **30 s**.

Spin eluates down and cover elution plate.

---

---

\* Reduction of atmospheric pressure

## 6 Appendix

### 6.1 Troubleshooting

Problem	Possible cause and suggestions
Low DNA yield	<p data-bbox="333 312 654 335"><i>Low DNA content of the sample</i></p> <ul data-bbox="333 357 983 639" style="list-style-type: none"> <li data-bbox="333 357 983 459">• The content of cell-free DNA in human plasma may vary over several orders of magnitude. DNA contents in the range of 0.1–1000 ng DNA per mL of plasma have been reported (see remarks in section 2.3).</li> <li data-bbox="333 481 983 639">• If the DNA concentration is measured with double strand specific dyes, e.g., PicoGreen®, make sure not to heat the DNA before measurement. Due to denaturation of DNA during the heat incubation step and the double strand specificity of certain DNA dyes, e.g., PicoGreen®, results may be inaccurate.</li> </ul>
Column clogging	<p data-bbox="333 662 777 684"><i>Sample contains residual cell debris or cells</i></p> <ul data-bbox="333 707 945 783" style="list-style-type: none"> <li data-bbox="333 707 945 783">• The plasma sample may have contained residual cells or cell debris. Make sure to use only clear plasma samples (see remarks in section 2.4).</li> </ul>
Discrepancy between $A_{260}$ quantification values and PCR quantification values	<p data-bbox="333 807 687 829"><i>Silica abrasion from the membrane</i></p> <ul data-bbox="333 852 972 1201" style="list-style-type: none"> <li data-bbox="333 852 972 1201">• Due to the typically low DNA content in plasma and the resulting low total amount of isolated DNA, DNA quantification via <math>A_{260}</math> absorption measurement is often hampered due to the low sensitivity of the absorption measurement. When performing absorption measurements close to the detection limit of the photometer, the measurement may be influenced by minor amounts of silica abrasion. In order to prevent incorrect <math>A_{260}</math>-quantification of small DNA amounts, centrifuge the eluate for 30 s at <math>&gt; 11.000 \times g</math> and take an aliquot for measurement without disturbing any sediment. Alternatively, use a silica abrasion insensitive DNA quantification method (e.g., PicoGreen® fluorescent dye).</li> </ul>
Unexpected $A_{260}/A_{280}$ ratio	<p data-bbox="333 1225 934 1248"><i>Measurement not in the range of photometer detection limit</i></p> <ul data-bbox="333 1270 969 1401" style="list-style-type: none"> <li data-bbox="333 1270 969 1401">• In order to obtain a significant <math>A_{260}/A_{280}</math> ratio, it is necessary that the initially measured <math>A_{260}</math> and <math>A_{280}</math> values are significantly above the detection limit of the photometer used. An <math>A_{280}</math> value close to the background noise of the photometer will cause unexpected <math>A_{260}/A_{280}</math> ratios.</li> </ul>



## 6.2 Ordering information

Product	REF	Pack of
NucleoSpin® 96 DNA Plasma	740873.1	1
	740873.4	4
NucleoSpin® 96 DNA Plasma Core Kit	740874.1	1
	740874.4	4
24-Square-well Block U-Bottom	740448.4	4
	740448.24	24
Square-well Block (96 wells)	740481	4
	740481.24	24
Round-well Block with Cap-strips	740475	4
	740475.24	24
Lysis Buffer PML	740835.125	125 mL
Binding Buffer PMB	740836.250	250 mL
Liquid Proteinase K	740396	5 mL
NucleoVac 96 Vacuum Manifold	740681	1
NucleoVac Vacuum Regulator	740641	1

## 6.3 Product use restriction/warranty

**NucleoSpin® 96 DNA Plasma** kit components are intended, developed, designed, and sold FOR RESEARCH PURPOSES ONLY, except, however, any other function of the product being expressly described in original MACHEREY-NAGEL product leaflets.

MACHEREY-NAGEL products are intended for GENERAL LABORATORY USE ONLY! MACHEREY-NAGEL products are suited for QUALIFIED PERSONNEL ONLY! MACHEREY-NAGEL products shall in any event only be used wearing adequate PROTECTIVE CLOTHING. For detailed information please refer to the respective Material Safety Data Sheet of the product! MACHEREY-NAGEL products shall exclusively be used in an ADEQUATE TEST ENVIRONMENT. MACHEREY-NAGEL does not assume any responsibility for damages due to improper application of our products in other fields of application. Application on the human body is STRICTLY FORBIDDEN. The respective user is liable for any and all damages resulting from such application.

DNA/RNA/PROTEIN purification products of MACHEREY-NAGEL are suitable for IN VITRO-USES ONLY!

ONLY MACHEREY-NAGEL products specially labeled as IVD are also suitable for IN VITRO-diagnostic use. Please pay attention to the package of the product. IN VITRO-diagnostic products are expressly marked as IVD on the packaging.

IF THERE IS NO IVD SIGN, THE PRODUCT SHALL NOT BE SUITABLE FOR IN VITRO-DIAGNOSTIC USE!

ALL OTHER PRODUCTS NOT LABELED AS IVD ARE NOT SUITED FOR ANY CLINICAL USE (INCLUDING, BUT NOT LIMITED TO DIAGNOSTIC, THERAPEUTIC AND/OR PROGNOSTIC USE).

No claim or representations is intended for its use to identify any specific organism or for clinical use (included, but not limited to diagnostic, prognostic, therapeutic, or blood banking). It is rather in the responsibility of the user or - in any case of resale of the products - in the responsibility of the reseller to inspect and assure the use of the DNA/RNA/protein purification products of MACHEREY-NAGEL for a well-defined and specific application.

MACHEREY-NAGEL shall only be responsible for the product specifications and the performance range of MN products according to the specifications of in-house quality control, product documentation and marketing material.

This MACHEREY-NAGEL product is shipped with documentation stating specifications and other technical information. MACHEREY-NAGEL warrants to meet the stated specifications. MACHEREY-NAGEL's sole obligation and the customer's sole remedy is limited to replacement of products free of charge in the event products fail to perform as warranted. Supplementary reference is made to the general business terms and conditions of MACHEREY-NAGEL, which are printed on the price list. Please contact us if you wish to get an extra copy.

There is no warranty for and MACHEREY-NAGEL is not liable for damages or defects arising in shipping and handling (transport insurance for customers excluded), or out of accident or improper or abnormal use of this product; defects in products or components not manufactured by MACHEREY-NAGEL, or damages resulting from such non-MACHEREY-NAGEL components or products.

MACHEREY-NAGEL makes no other warranty of any kind whatsoever, and SPECIFICALLY DISCLAIMS AND EXCLUDES ALL OTHER WARRANTIES OF ANY KIND OR NATURE WHATSOEVER, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, INCLUDING, WITHOUT LIMITATION, AS TO THE SUITABILITY, REPRODUCTIVITY, DURABILITY, FITNESS FOR A PARTICULAR PURPOSE OR USE, MERCHANTABILITY, CONDITION, OR ANY OTHER MATTER WITH RESPECT TO MACHEREY-NAGEL PRODUCTS.

In no event shall MACHEREY-NAGEL be liable for claims for any other damages, whether direct, indirect, incidental, compensatory, foreseeable, consequential, or special (including but not limited to loss of use, revenue or profit), whether based upon warranty, contract, tort (including negligence) or strict liability arising in connection with the sale or the failure of MACHEREY-NAGEL products to perform in accordance with the stated specifications. This warranty is exclusive and MACHEREY-NAGEL makes no other warranty expressed or implied.

The warranty provided herein and the data, specifications and descriptions of this MACHEREY-NAGEL product appearing in MACHEREY-NAGEL published catalogues and product literature are MACHEREY-NAGEL's sole representations concerning the product and warranty. No other statements or representations, written or oral, by

MACHEREY-NAGEL's employees, agent or representatives, except written statements signed by a duly authorized officer of MACHEREY-NAGEL are authorized; they should not be relied upon by the customer and are not a part of the contract of sale or of this warranty.

Product claims are subject to change. Therefore please contact our Technical Service Team for the most up-to-date information on MACHEREY-NAGEL products. You may also contact your local distributor for general scientific information. Applications mentioned in MACHEREY-NAGEL literature are provided for informational purposes only. MACHEREY-NAGEL does not warrant that all applications have been tested in MACHEREY-NAGEL laboratories using MACHEREY-NAGEL products. MACHEREY-NAGEL does not warrant the correctness of any of those applications.

Last updated: 07/2010, Rev. 03


Please contact:  
MACHEREY-NAGEL GmbH & Co. KG  
Tel.: +49 24 21 969-270  
tech-bio@mn-net.com

---

**Trademarks/Disclaimer:**

Cell-Free DNA BCT is a registered trademark of Streck, Inc.  
NucleoSpin is a registered trademark of MACHEREY-NAGEL GmbH & Co. KG  
Picogreen is a registered trademark of Molecular Probes, Inc.

All used names and denotations can be brands, trademarks, or registered labels of their respective owner – also if they are not special denotation. To mention products and brands is only a kind of information (i.e., it does not offend against trademarks and brands and can not be seen as a kind of recommendation or assessment). Regarding these products or services we can not grant any guarantees regarding selection, efficiency, or operation.



[www.mn-net.com](http://www.mn-net.com)

## MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0  
 Fax: +49 24 21 969-199  
 E-mail: [info@mn-net.com](mailto:info@mn-net.com)

CH:

Tel.: +41 62 388 55 00  
 Fax: +41 62 388 55 05  
 E-mail: [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)

FR:

Tel.: +33 388 68 22 68  
 Fax: +33 388 51 76 88  
 E-mail: [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)

US:

Tel.: +1 484 821 0984  
 Fax: +1 484 821 1272  
 E-mail: [sales-us@mn-net.com](mailto:sales-us@mn-net.com)