

MACHEREY-NAGEL

Plasmid DNA

desalting and concentration

Bioanalysis



### NucleoSnap Finisher Midi / Maxi

Speed up your NucleoBond® plasmid preparations

- No time-consuming isopropanol precipitation
- Process 12 samples in less than 10 minutes
- 100% plasmid DNA recovery



**NEW**  
column design  
for VACUUM  
processing

**MACHEREY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)



Since 1911

# NucleoSnap Finisher Midi / Maxi

## The fastest way to desalt and concentrate DNA after your NucleoBond® preparation

Eluates from anion-exchange plasmid preparations (such as NucleoBond® preparations) contain high salt concentrations. MACHEREY-NAGEL has developed a new vacuum-driven accessory that enables desalting and concentrating plasmid DNA in a few minutes without any loss of yield or quality.

No time-consuming centrifugation steps necessary!



### Product at a glance

Technology	Silica-membrane technology
Format	Snap-off column, vacuum processing
Compatibility	Eluates from NucleoBond® preparations (Xtra Midi and Maxi, PC 20, PC 100, PC 500)
Sample material	DNA eluates from e.g., anion-exchange purification kits
Typical yield	Dependent on starting material
Elution volume	≥100 µL
Preparation time	<10 min/12 preps
Binding capacity	1.5 mg

Also available for centrifugation: NucleoSpin® Finisher Midi (see ordering information)



### Application data



#### Ultra-fast plasmid DNA desalting and concentration

Eluates from NucleoBond® Xtra Maxi preparations were precipitated with the NucleoSnap Finisher Maxi and a conventional isopropanol precipitation. 6 preparations can be processed in only 7 minutes with the NucleoSnap Finisher Maxi.

#### 100% DNA recovery

DNA yield in eluates from a NucleoBond® Xtra Maxi preparation were measured before and after precipitation with the NucleoSnap Finisher Maxi and isopropanol. With the NucleoSnap Finisher, 100% DNA recovery can be achieved.

### Ordering information

Product	Specifications	Preps / pack of	REF
NucleoSnap Finisher Midi	Vacuum-processed precipitation of eluates from NucleoBond® PC 20, PC 100 and Xtra Midi (EF) preparations	10/50	740434.10/50
NucleoSnap Finisher Maxi	Vacuum-processed precipitation of eluates from NucleoBond® PC 500 (EF) and Xtra Maxi (EF) preparations	10/50	740435.10/50
NucleoVac 24 Vacuum Manifold	Vacuum manifold for processing 1–24 NucleoSpin® or NucleoSnap columns - manifold, NucleoVac Mini Adapter, Luer plugs, tubing connection, closing plug	1	740299
NucleoSpin® Finisher Midi	Precipitation of plasmid DNA eluates from NucleoBond® PC 100 and Xtra Midi (EF) preparations, processing by centrifugation	10/50	740439.10/50

NucleoBond® and NucleoSpin® are registered trademarks of MACHEREY-NAGEL GmbH & Co. KG · Düren · Germany

www.mn-net.com

# MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

DE / International:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

CH:

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

FR:

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

US:

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: sales-us@mn-net.com

