

Dye Terminator Removal: NucleoSEQ

Unincorporated dye terminators will negatively affect analysis of sequencing results. Excess of dye terminators causes so-called “dye blobs” resulting in a partly unreadable sequence. **NucleoSEQ** will remove unincorporated dye terminators. The subsequent analysis is of high quality with long reading length and minimized background.

Your reasons to use NucleoSEQ

- ⇒ gel-filtration technology optimized for efficient removal of dye terminators, e.g. BigDye™ Terminators
- ⇒ convenient single spin columns
- ⇒ time saving, no ethanol precipitation necessary
- ⇒ long-term storage at room temperature
- ⇒ cost efficient alternative to competitive products

Cleanup of sequencing reactions with NucleoSEQ columns ensures high-quality sequencing results.



Sequencing profile of plasmid DNA (pGEM®-T Easy). Plasmid DNA was purified using **NucleoSpin® Plasmid**. Sequencing reaction was performed with ABI PRISM® BigDye™ Terminator Cycle Sequencing kit, purified with **NucleoSEQ**, and analyzed on an ABI 310 sequencer.

BIOKÉ
sharing knowledge

Plesmanlaan 1d
2333 BZ Leiden
The Netherlands
T. +31 (0)71 568 10 00
T. Belgium: 0800 71640
F. +31 (0)71 568 10 10
info@bioke.com
www.bioke.com

Principle

NucleoSEQ columns are designed for fast, effective and cost efficient clean-up of sequencing reactions. The spin columns are prefilled with a dry size exclusion matrix which allows an efficient removal of dye terminators, e.g. BigDye™ Terminators: The gel-filtration material consists of spheres with uniform pores and separates molecules according to molecular weight. After applying the sequencing reaction to the **NucleoSEQ** column, small dye terminators and other impurities e.g. salts, nucleotides, primers, traces of organic solvents are retained into the pores while labeled DNA fragments are excluded and recovered in the flow-through with high yield.

Handling

In order to achieve long-time storage life at room temperature, **NucleoSEQ** columns are prefilled with dry gel-filtration resin. The matrix can easily be hydrated by adding water followed by an incubation period (>30 min). Hydrated columns are ready to use and can be stored at 4°C for 14 days.

A first short centrifugation step removes remaining storage buffer. After loading the sample onto the column and a second centrifugation step, the DNA fragments of interest are recovered in the flow-through.

Receiver Columns 20 µm

Receiver columns are micro spin-columns with an inserted hydrophobic frit of 20 µm pore size. They can be used for general filtration purposes as well as for retaining chromatographic resins (e.g. **NucleoSil**® C18, Sephadex® G25, G50, or Sephacryl® S200). Receiver columns 20 µm are delivered with a closed outlet inserted into a collection tube and are for use with suitable bench-top centrifuges.

- ⇒ Filtration of viscous solutions
- ⇒ Filtration of swabs, e.g. buccal swabs
- ⇒ Desalting of protein solutions
- ⇒ ...

Trademarks:


NucleoSpin® and **NucleoSil**® are registered trademarks of MACHERY-NAGEL GmbH & Co. KG (Germany). pGEM® is a registered trademark of Promega Corporation (USA). ABI PRISM® is a registered trademark of Applied Biosystems (USA). Sephadex® and Sephacryl® are registered trademarks of Pharmacia AB, Uppsala (Sweden). BigDye™ is a trademark of Applied Biosystems (USA).

Ordering Information:


| Product | Cat. No. |
|--|------------------------|
| NucleoSEQ (10/50/250 preps) | 740523.10 / .50 / .250 |
| Receiver Columns 20 µm (10/50/250 columns) | 740522.10 / .50 / .250 |

Distributed by:

For more information regarding the use of MN products, please contact your local representative or visit MN directly under www.mn-net.com.



NucleoSEQ procedure



NUSQO_001/15.019.2005 PD
Printed in Germany

MACHERY-NAGEL

MACHERY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8 · D-52355 Düren · Germany

Germany

Tel.: +49 (0) 2421 969-275
e-mail: tech-bio@mn-net.com

USA

Tel.: +1 610-559-9848
e-mail: sales-us@mn-net.com

France

Tel.: +33 (0) 388 682268
e-mail: sales-fr@mn-net.com

Switzerland

Tel.: +41 (0) 62 388 55 00
e-mail: sales-ch@mn-net.com

