

Plate Compatibility Table

Plate compatibility information is for guidance only. Please check installed heat block.

Please refer to instruments manual for information or contact our customer support.

For compatibility information of the Vari-Plates™ refer to the table entry for respective frame.

Framestar® 2 component

Standard



FrameStrip™
4ti-0785/0786



Break-A-Way
4ti-1000



Low Profile Break-A-Way
4ti-1200



384 well 480
4ti-0380/0381



384 well
4ti-0384



96 well 480
4ti-0950/0951



96 well FastPlate
4ti-0910/C



96 well semi-skirted
(cut H1) 4ti-0900



96 well semi-skirted
(cut A12) 4ti-0770



96 well skirted
4ti-0960



96 well with upstand
4ti-0730



96 well non-skirted,
low profile 4ti-0720



96 well non-skirted,
std profile 4ti-0710



96 well non-skirted
/ Tear-A-Way
4ti-0750 / 4ti-0750/TA



96 well semi-skirted
4ti-0760



96 well semi-skirted 480
4ti-0955



96 well skirted
4ti-0740



8-tube strip
4ti-0780/0784



Vari-Strip™
4ti-0753

	7000, 7300, 7500, 7700, 7900	7900 HT	FAST 7500, 7900, 7900HT	ABI Life Technologies	StepOne Plus™	QuantStudio™	ViiA7™	CFX 384	CFX96	BIO-RAD	iCycler™ IQ, IQ™ 4, IQ™ 5	MyiQ™	Chromo4™	BIO-RAD MJ RESEARCH	Opticon™, Opticon2™	Mastercycler™ ep realplex	EPENDORF	LC480	ROCHE	LC96	Nano	Mx3000PT™	STRATAGENE / AGILENT	Mx3005PT™	Mx4000™	TECHNE	Quantica
FrameStrip™ 4ti-0785/0786	•										•	•				•							•	•	•	•	•
Break-A-Way 4ti-1000	•										•	•				•							•	•	•	•	•
Low Profile Break-A-Way 4ti-1200									•					•	•	•											
384 well 480 4ti-0380/0381																		•									
384 well 4ti-0384	•					•	•	•																			
96 well 480 4ti-0950/0951																•											
96 well FastPlate 4ti-0910/C			•	•	•	•																					
96 well semi-skirted (cut H1) 4ti-0900									•	•				•								•					
96 well semi-skirted (cut A12) 4ti-0770	•					•	•				•	•				•						•	•	•	•	•	•
96 well skirted 4ti-0960									•				•	•	•												
96 well with upstand 4ti-0730	•	•				•	•									•									•		
96 well non-skirted, low profile 4ti-0720									•					•	•	•										•	
96 well non-skirted, std profile 4ti-0710											•	•				•						•	•	•	•		
96 well non-skirted / Tear-A-Way 4ti-0750 / 4ti-0750/TA											•	•		•		•						•	•	•	•		
96 well semi-skirted 4ti-0760	•					•	•									•								•	•	•	•
96 well semi-skirted 480 4ti-0955																		•									
96 well skirted 4ti-0740									•				•	•	•								•	•	•		
8-tube strip 4ti-0780/0784	•										•	•				•						•	•	•	•	•	•
Vari-Strip™ 4ti-0753																			•								

Plate Compatibility Table

Framestar® 2 component

Standard

	GeneAmp® 2400	GeneAmp® 2700, 2720, 9600	GeneAmp® 9700	ABI / LIFE TECHNOLOGIES GeneAmp® 9800 FAST Block Veriti (384-well block)	Veriti 0.1ml (96-well block) FAST Veriti 0.2ml (96-well block)	T1 Thermocycler	T3 Thermocycler	TGradient	BIOMETRA TPersonal Cycler TRobot, Uno Uno II	C1000	Genecycler	iCycler™	BIO-RAD MyCycler™ S1000 T100	Mini Gradient	Personal	PTC100™ (96-well block only)	DNA Engine™, DNA Dyad™, DNA Tetrad™	CORBETT RESEARCH (Qiagen) Palm Cycler (Qiagen) Palm Cycler 384	M384	EPPENDORF MasterCycler® MasterCycler® EP Gradient MasterCycler® Gradient	SWIFT Swift	ESCO
FrameStrip™ 4ti-0785/0786	•	•	•		•	•		•		•	•	•	•			•	•			•	•	
Break-A-Way 4ti-1000		•	•		•	•		•		•	•	•	•			•	•			•	•	
Low Profile Break-A-Way 4ti-1200						•		•	•	•	•	•	•			•	•	•				
384 well 480 4ti-0380																						
384 well 4ti-0384			•	•		•				•	•	•	•				•	•				
96 well 480 4ti-0950																						
96 well FastPlate 4ti-0910/C				•	•																	
96 well semi-skirted (cut H1) 4ti-0900						•		•	•	•	•	•	•		•	•	•				•	•
96 well semi-skirted (cut A12) 4ti-0770	•	•			•	•		•	•	•	•	•	•		•	•	•			•	•	
96 well skirted 4ti-0960 / 4ti-0757						•		•	•	•	•	•	•		•	•	•				•	•
96 well with upstand 4ti-0730	•	•			•																	
96 well non-skirted, low profile 4ti-0720						•		•	•	•	•	•	•		•	•	•			•	•	
96 well non-skirted, std profile 4ti-0710	•	•				•		•	•	•	•	•	•		•	•	•			•	•	•
96 well non-skirted / Tear-A-Way 4ti-0750 / 4ti-0750/TA	•	•				•		•	•	•	•	•	•		•	•	•			•	•	
96 well semi-skirted 4ti-0760	•	•			•	•		•	•	•	•	•	•		•	•	•			•	•	
96 well semi-skirted 480 4ti-0955																						
96 well skirted 4ti-0740						•		•	•	•	•	•	•		•	•	•				•	•
8-tube strip 4ti-0780/0784	•					•	•	•	•	•	•	•	•		•	•	•			•	•	•
Vari-Strip™ 4ti-0753			•	•	•	•		•	•	•	•	•	•		•	•	•			•	•	•