

Your personal single cell sorters are here!



Namo™
Single Cell Dispenser
up to 3 colors



Hana™
Single Cell Dispenser
scatter + 2 colors



**Fast single cell
dispensing**
in 1 min/96-well plate



Bulk sorting
sorts 100M
cells in 5 min



Sterile sorting
with disposable
cell cartridge



**Easy with
automated**
initialization in 2 min



Low pressure
to preserve
cell viability


Namocell
The Fastest Way to a Single Cell

BIOKÉ
www.bioke.com